



HINT 365

HealthyInstantTracker

User guide

The blood-pressure diary: record readings by voice or with a photo, and prepare the report for the doctor.

1. What HINT 365 is

HINT 365 is a **diary**: it records the values of your blood-pressure monitor, **SYS** (systolic), **DIA** (diastolic) and **PUL** (pulse), with date and time, and shows them in tables, charts and reports.

HINT 365 is not a medical device and gives no medical advice. It does not assess the values, does not say whether they are normal or high, makes no diagnosis and does not replace the doctor. The assessment is always the doctor's.

2. Free, and optional: what is paid

Feature	Cost
Recording readings by voice, tables, charts, PDF and Excel reports, reminders, recovery on a new phone	Free
Scanning: reading the numbers from a photo of the monitor with artificial intelligence	Optional, paid

Scanning uses **Anthropic's** artificial intelligence and costs about **half a US cent per photo**. You pay it, with your own prepaid Anthropic credit: nobody else pays for your photos and you do not pay for theirs. Without it, everything else in the app works.

3. Installing

1. On the Android phone open <https://personalhealthy-api.comevava76.workers.dev/download> and download **HINT.apk**.
2. Open the file: Android asks to allow installs from this source. Allow it.
3. Open HINT 365, read the note about your data, tick "**I have read and accept**" and tap "**Sign in with Google**".
4. From then on the app opens with the phone's **fingerprint, face or screen lock**.

Google is used only to create your diary and to find it again on a new phone. It is not asked day to day.

4. Recording a reading

Record (free)

1. Tap **Record** and say the three values in order: **SYS**, **DIA**, **PUL** (e.g. "127, 80, 70"). Pauses are fine.
2. The app stops by itself once it has heard three numbers and shows them.
3. Check them and tap **Save**. Date and time are those of the moment you spoke.

The app refuses impossible values (for example DIA higher than SYS, a missing number) and marks unusual ones in amber, to check again.

Scan (optional, paid)

1. Tap **Scan** and photograph the monitor's display: close, straight, without reflections.
2. The app reads the numbers and shows them. Check they match the display and tap **Save**.

5. Turning on scanning (AI credit)

At first start the app shows a welcome screen explaining the two ways: you can tap "**Continue without scanning**" and turn it on later from **Admin** → **Token** → "**Add my key**". There are four steps, each with its button:

1. **Create an Anthropic account** (free, with your email).
2. **Add credit**: load a small prepaid credit (a few dollars last for months).
3. **Create a key**: it starts with **sk-ant-**; copy it.
4. **Save the key here**: back in the app, tap "**Paste the key**", type the credit you see on Anthropic and tap **Save**.

The key is tested and stored encrypted on the server; it is never shown again, not even to you. You can replace or delete it at any time. The balance in the app is an estimate: the real one is on Anthropic.

6. The screens

Tab	What you find
Blood pressure	The Record and Scan buttons, the last reading, the last 7 days table (one average per day) and the summary with highest values and averages.
Report	7, 15 or 30 days: daily chart, period averages, Day average and Evening average , extreme values, sending the PDF for the doctor and the Excel file.
Admin	Credits (estimated balance, recharge), Token (your key), Readings database (delete one reading or all), Identity (Google account, share the app, delete the account).

7. Reminders

Every day at **9:00** and **17:00** "Remember to measure your blood pressure" arrives, with a number on the icon. At midnight they are cleared. Turn them off in the phone settings ("Measurement reminders").

8. New or lost phone

Install HINT 365 on the new phone and tap "**Sign in with Google**" with the same account: all your readings are there. The old phone is disconnected. To update the app, install the new version **over** the existing one, without uninstalling.

9. Your data

- **What is saved:** the email and identifier of your Google account (only to find your diary again), your readings with date and time and their source (photo or voice), and your encrypted Anthropic key if you add it.
- **Photos and voice:** the photo is sent to Anthropic only to read the numbers and is not kept by the app. Speech is transcribed by the phone; only the three numbers reach the app.
- **Where:** on **Cloudflare** servers, with the database in the **European Union**. The app's code is on **GitHub**.
- **Who sees it:** only you. Other users do not. The person who runs the service can reach the database only for maintenance.
- **Deleting:** in Admin you can delete one reading, all readings, or the account with all its data. It is immediate and final.

Full privacy policy: <https://personalhealthy-api.comevava76.workers.dev/privacy>

10. Responsibility

No medical advice. HINT 365 collects and shows the data you enter. Every health decision is to be made with the doctor. The operator of HINT 365 is not responsible for medical choices made on the basis of the app's data.

Always check the values. Photo reading and speech recognition can be wrong: before saving, compare the numbers with the monitor's display.

Service without warranties. HINT 365 is offered free of charge, as it is, with no guarantee of continuous availability. Data is kept on third-party cloud infrastructure (Cloudflare) and the code on GitHub. To the extent permitted by law, the operator is not liable for service interruptions, errors or data loss.

Keep a copy. Regularly export the PDF or Excel file from the Report tab and keep it.

This document describes the service and is not legal advice.