



HINT 365

HealthyInstantTracker

User guide

The blood-pressure diary: record readings by voice or with a photo, and prepare the report for the doctor.

1. What HINT 365 is

HINT 365 is a **diary**: it records the values of your blood-pressure monitor, **SYS** (systolic), **DIA** (diastolic) and **PUL** (pulse), with date and time, and shows them in tables, charts and reports.

HINT 365 is not a medical device and gives no medical advice. It does not assess the values, never says whether they are normal or high, makes no diagnosis and does not replace your doctor. Assessment is always up to the doctor.

2. Costs, plainly

What	Cost	Paid to
HINT 365: the app, charts, reports, Web Dashboard, keeping your readings	Free today. Once it is on the Play Store: 5 \$ a year, cancel any time	The owner of the app, through Google Play
AI features, optional: Scan of the monitor (available now); local lab report import (no Anthropic fee)	Pay-per-use, about 0.006 \$ a photo	Anthropic, from your own credit

Without the AI features you pay nothing more. The Anthropic credit is yours: you manage it on Anthropic, it goes down only when you use the AI, and you can use it for other things too. The owner of HINT 365 receives none of it.

3. Installing and updating

1. On your Android phone open <https://personalhealthy-api.comevava76.workers.dev/download> and download **HINT.apk**.
2. Open the file: Android asks to allow installing from this source. Allow it.
3. Tap "**Sign in with Google**": sign-in uses the OAuth 2.0 and OpenID Connect standards, and HINT 365 never receives your password.
4. Read the **terms of use** and accept them: they appear at the first installation and after **every update**. Without acceptance the app does not open.

5. From then on the app opens with the phone's **fingerprint, face or screen lock**. If the phone has no screen lock at all, the app opens unprotected and Blood pressure shows an orange notice that leads to the settings: protecting the phone is your responsibility.

To update, install the new version **over** the existing one, without uninstalling. If an old version is switched off, the app shows only the button to download the latest: your data stay.

4. Recording a reading

Record

1. Tap **Record** and say the three values in order: **SYS**, **DIA**, **PUL** (e.g. "127, 80, 70"). You can pause.
2. The app stops by itself once it has heard three numbers and shows you the values.
3. Check them and tap **Save reading**. Date and time are those of the moment you spoke.

The app refuses impossible values (for example DIA higher than SYS, a missing number) and marks unusual ones in orange, to check again.

✦ **Scan (AI feature)**

1. Tap **Scan** and photograph the monitor's display: close, straight, without glare.
2. The app reads the numbers and shows them. Check they match the display and tap **Save reading**. The photo is not kept.

Scan is on only when Anthropic confirms that the key is valid and there is credit. The check is automatic: when the app opens, after an error and when you pull the screen down. If the credit runs out, Scan switches off and the app tells you.

5. Turning on the AI features

Tap ✦ **Turn on AI** (at the top of the Blood pressure screen, or in Admin). There are four steps, each with its own button:

1. **Create an Anthropic account** (free, with your email).
2. **Add credit**: load a small prepaid credit (a few dollars last for many photos).
3. **Create a key**: it starts with **sk-ant-**; copy it.
4. **Save the key here**: back in the app, paste it and tap **Save**.

The key is tested with Anthropic and stored encrypted on the server; it is never shown again, not even to you. You can replace or delete it at any time. Anthropic lets no app read the balance: the app shows your **estimated spending**, and the exact balance is on your Anthropic page.

6. The screens

At the bottom every tab has an icon and its name. Only the owner of the app also sees **Owner** (shield ↗), which opens the owner's area on the web.

Tab	What you find
Blood pressure (heart)	The Record and Scan buttons, the last reading, the week (one mean per day) and the summary; at the bottom Send PDF , Send Excel and Web Dashboard .
Lab results (sheet with a flask)	Import lab report , the progress bar while it uploads and below it the upload history , day by day, with the outcome of each (saved, already imported, not saved and why). The results are in the Web Dashboard.
Dashboard (screen with a chart ↗)	Not a screen of the app: it opens the Web Dashboard in the browser, on the Blood pressure tab. The last 7 days: chart with one dot per day (the mean of the day), means, morning and evening with the balance, extremes, pulse; in the tables the moment is AM (before 12:00) or PM . Then the Lab results tab.
Admin	Credits (estimated spending, link to the balance on Anthropic, Recharge), AI features (state, replace or delete the key), Costs , Readings database (delete one reading or all), Identity (Google account, sign out, delete the account) and, at the bottom, Share the app .

7. Web Dashboard and the doctor

1. Tap **Dashboard** in the bar at the bottom (or **Web Dashboard** at the end of Blood pressure and Lab results): the browser opens already signed in to your account (the link lasts 60 seconds).
2. The first time, accept the **technical cookie** that keeps you signed in for 7 days. It is the only cookie: no advertising, statistics or third parties. The **Cookies** link at the bottom shows your choice and changes it; there is no Sign out button: the session ends by itself.
3. You see large charts of the last 7 days and the list of readings, and you can download the same A4 PDF as the app.
4. **Send to the doctor**: choose **Email** or **WhatsApp**. The PDF of the report and a read-only link valid 7 days go together, with a ready message in the phone's language ("Good morning, I am sending the blood pressure report recorded in the period ..."; email subject "Blood pressure report"). On a phone the PDF is attached for you; on a computer it is downloaded and you attach it. You can withdraw the links at any time.

8. Reminders

Every day at **9:00** and **17:00** "Remember to measure your blood pressure" arrives, with a number on the icon. They reset at midnight. Turn them off in the phone's settings ("Measurement reminders").

9. New or lost phone

Install HINT 365 on the new phone and tap **"Sign in with Google"** with the same account: you find all your readings. The old phone is signed out.

10. Your data

- **What is saved**: your readings with date, time and source (voice or photo); of your Google account only an encrypted fingerprint, **never the email**, which stays on your phone; your Anthropic key,

encrypted, if you add it.

- **Photo and voice:** the photo goes to Anthropic only to read the numbers and is not kept. Speech is transcribed by the phone; only the three numbers reach the server.
- **Where and for how long:** on **Cloudflare**, database in the **European Union**. Readings are kept for at most **365 days**: older ones are deleted every night. Export PDF or Excel to keep them longer.
- **Who sees them:** only you, and the doctor with the link you send. The owner of the app sees only usage numbers tied to an anonymous code (when you joined, last opened, how many readings), never values or reports. Technical errors go into a short error log, without reading values, kept 90 days, to fix problems.
- **Deletion:** in Admin you can delete one reading, all readings or the account. It is immediate and final; only the proof of the accepted terms remains.

Full privacy policy: <https://personalhealthy-api.comevava76.workers.dev/privacy> · Terms of use: <https://personalhealthy-api.comevava76.workers.dev/terms>

11. Responsibility

No medical advice. HINT 365 collects and shows the data you enter. Every health decision is to be taken with your doctor.

Always check the values. Photo reading and speech recognition can be wrong: before saving, compare the numbers with the monitor's display.

No warranty. HINT 365 is offered as is, without a guarantee of continuous availability; the terms of use accepted in the app apply.

Keep a copy. Regularly export the PDF or the Excel file (at the bottom of the Blood pressure tab) and keep it.

This document describes the service and is not legal advice.

Lab results

Only your own reports. One account is one person: upload only your own reports and readings, never those of a relative. HINT 365 does not know whose data they are and cannot separate them (terms of use).

1. Open **Lab results** → **Import lab report** and choose the laboratory's PDF or a photo of the report (20 MB maximum; unencrypted PDFs up to 15 pages). Other kinds of file cannot be chosen.
2. The app does the rest, also if you switch tab: a bar fills through **Processing** → **Scanning** → **Uploading**; the outcome stays a few seconds, then the bar goes away and the outcome is in the upload history. The document is read only on the phone; the original, patient identifiers and full text are not uploaded.
3. It is saved only when everything is clear: one report date and every result line understood, also tests outside our catalog (with the name as printed). If anything is unclear nothing is saved and the app says why.

4. Never twice: the same file, or the same test on the same date, is not saved again. An older report is saved with its own date. Results are never deleted by themselves: they stay until you delete them. The document stays on the phone; only the date and the values as printed reach the server, encrypted in transit, never a name or any reference to the person.
5. On the web, **Lab results** shows a table: one row per test, one column per report date; a dash where a test is missing. In orange with ↑ or ↓: a value outside the reference printed on the report (a comparison only, not an assessment). **Delete** under a date removes that column; the report can then be imported again. **Delete all lab results**, under the table, removes them all. **Share** and **PDF** give the same table.

The app and web follow the device/browser language: Italian, English, German or French (English fallback). Lab import needs no Anthropic key and has no Anthropic charge. Medical interpretation belongs to your doctor.